

INVERCLYDE BIOLOGICALS

Tel: 01698 749368 Fax: 01698 740521

COSHH SUBSTANCES INFORMATION SHEET FOR A FINISHED PRODUCT

SUBSTANCE	LIQ-A-TEC KIT
	Alternative name: Amniotic Fluid Detection Kit
Exposure limits;	Long term (8hr.) N/A ppm mg./cu.m. Short term (10min.) N/A ppm mg./cu.m.
=====	
HAZARD;	Only minor flammable problem with the reagent
=====	
FIRST AID;-	
Eyes;-	Irrigate thoroughly with water for 10 min. If discomfort persists SEEK MEDICAL ADVICE.
Lungs;-	Remove from exposure, rest and keep warm. SEEK MEDICAL ADVICE.
Skin;-	Drench with plenty of water. Remove contaminated clothing.
Mouth;-	Wash out thoroughly with water. Give plenty of water as a mouth wash. Mouth wash water NOT to be swallowed. DO NOT INDUCE VOMITING SEEK MEDICAL ADVICE if the liquid has been swallowed.
=====	
FIRE FIGHTING;-	N/A
=====	
SPILLAGE;-	Mop up with plenty of water and run to waste diluting greatly with running water.
=====	
PRECAUTIONS DURING USE;-	
S24/25	Avoid contact with skin and eyes.
S56	Mouth pipetting is expressly forbidden.
=====	
PROTECTIVE MEASURES;-	
Respirator	N/A
Ventilation	N/A
Gloves	Rubber or latex.
Eye protection	Not necessary
Other Measures	Always wash hands after using this reagent
=====	
STORAGE AND HANDLING;-	
Special requirements	None
=====	
CONCLUSION;-	THIS REAGENT PRESENTS ONLY A MINOR FLAMMABLE RISK TO STAFF
=====	

If further information is required please contact INVERCLYDE *Biologicals*, Bellshill, Scotland.
Telephone: 01698 749368 Fax: 01698 740521